

# Applications of Nanorobots in Medical Technique

**Bhushan Bhakare**

Department of Pharmacy

Samarth College of Pharmacy, Belhe, Pune, Maharashtra, India

bhushanbhakare0843@gmail.com

**Abstract:** *Nanorobots, with their minuscule size ranging from 1 to 100 nanometers, have emerged as a revolutionary force in the realm of medical technology. This abstract explores the diverse applications of nanorobots, showcasing their potential to reshape diagnostics, drug delivery, surgery, and patient monitoring. Equipped with advanced sensors and imaging capabilities, nanorobots hold the promise of early disease detection by navigating the bloodstream and monitoring biomarkers. Their role in drug delivery involves precise transport of therapeutic agents to targeted cells, minimizing side effects and enhancing treatment efficacy. In the surgical domain, nanorobots enable minimally invasive procedures, performing intricate tasks at the cellular or molecular level. Continuous monitoring and intervention capabilities make nanorobots invaluable for proactive healthcare management. Moreover, in targeted cancer therapy, these microscopic marvels can deliver therapeutic payloads directly to cancer cells, offering a promising avenue for minimizing collateral damage. Despite these transformative potentials, challenges related to biocompatibility and safety necessitate further exploration. The applications of nanorobots in medical techniques represent a paradigm shift, holding the potential to redefine patient care and treatment strategies.*

**Keywords:** Nanorobots, Nanotechnology, Cancer Treatment, Surgery, Medicine, Healthcare, Gene Therapy, Component, Nano-surgery, Disease Management

## I. INTRODUCTION

Nanotechnology, an interdisciplinary field encompassing chemistry, physics, materials science, and biology, integrates diverse expertise to advance novel technologies. Within this expansive realm, nanorobots emerge as noteworthy nano-devices designed to offer protection or treatment against pathogens in humans<sup>1</sup>. Operating at nanoscale dimensions ranging from 1 to 100 nanometers, nanorobots exhibit precision in executing specific tasks. These microscopic entities hold the potential to operate at atomic, molecular, and cellular levels, contributing to transformative applications in both medical and industrial domains. In medical contexts, nanorobots are envisioned to play a pivotal role in targeted protection and treatment, showcasing their versatility in addressing challenges at the forefront of healthcare and industrial processes<sup>2</sup>. medicine, bioinformatics, and computers can lead to the development of the nanorobot drug delivery system. Some of the examples are respirecyte nanorobots, microbivore nanorobots, surgical nanorobots, and cellular repair nanorobots<sup>3</sup>. They have a diameter of about 0.5 to 3 microns and will be constructed out of parts with dimensions in the range of 1 to 100 nanometers<sup>4</sup>. The main element used in nanorobots is carbon because of its inertness and strength in the form of diamond and fullerene. They generally have an exterior passive diamond coating to avoid attack by the host immune system<sup>3</sup>. Techniques like Scanning Electron Microscopy (SEM) and Atomic Force Microscopy (AFM) are being employed to understand the molecular structure of the nanoscaled devices<sup>5,6</sup>.

### Advantages of Nanorobots over Conventional Medical Techniques:

- Precision Targeting: Nanorobots operate at molecular and cellular levels, providing precise targeting for diagnostics, drug delivery, and surgery<sup>7</sup>.
- Minimally Invasive: Nanorobots enable minimally invasive procedures, reducing recovery times and lowering the risk of complications compared to conventional surgeries<sup>7</sup>.



- **Early Disease Detection:** Equipped with advanced sensors, nanorobots can detect diseases at an early stage, allowing for prompt intervention and personalized treatment.
- **Targeted Drug Delivery:** Nanorobots deliver therapeutic agents directly to specific cells, minimizing side effects associated with systemic drug administration<sup>7,8</sup>.
- **Continuous Monitoring:** Nanorobots offer continuous monitoring of physiological parameters, enabling proactive healthcare management and timely intervention.
- **Customization:** Nanorobots facilitate personalized medicine by tailoring treatments to individual patients based on their unique biological characteristics.
- **Enhanced Imaging:** Equipped with advanced imaging technologies, nanorobots provide detailed and real-time information for improved diagnostics.
- **Reduced Side Effects:** Precise targeting in drug delivery minimizes side effects, enhancing the overall safety and efficacy of treatments.

Scientists and researchers are working on a more robust, reliable, and bio-compatible approach. Instead of curing of the outside, they plan to defend the body from the inside. That is where medical nanorobots come in. The major advantages of this technology provide are

1. Minimal or no tissue trauma.
2. Considerably less recovery time.
3. Less post-treatment care required.
4. Continuous monitoring and diagnosis from the inside.
5. Rapid response to a sudden change.

Some added features of nanorobotics. Would also enable us to do the following:

1. It can store and process previous data, identify patterns, and hence, help to predict the onset of an ailment.
2. It can guide externally or as per program, targeting specific locations.
3. It has the ability to deliver payloads such as drugs, or healthy cells to the specific site<sup>9,10</sup>.

#### **Disadvantages of nanorobots**

- **Biocompatibility Risks:** Concerns about the compatibility of nanorobots with biological systems, leading to potential immune responses or toxicity.
- **Ethical Concerns:** Ethical dilemmas surrounding privacy, autonomy, and misuse of nanorobot technology.
- **Safety Issues:** Risks of malfunctions, unintended side effects, or the release of toxic materials in the human body.
- **Complex Manufacturing:** Challenges in ensuring consistent and precise manufacturing of nanorobots at scale.
- **Limited Power Supply:** Constraints on sustaining reliable power sources for nanorobots, especially within the human body.
- **Navigational Challenges:** Difficulty in navigating nanorobots to specific target locations within complex biological environments.
- **Regulatory Hurdles:** Evolving regulatory frameworks and standards for ensuring the safety and efficacy of nanorobotic applications.
- **Cost Considerations:** High expenses associated with research, development, and integration of nanorobotic technologies<sup>11,12</sup>.

#### **Different Approaches to Nanorobotics:**

Nanorobots, also known as nanobots, can be designed and implemented using various approaches, each tailored to specific applications and functionalities. Some notable approaches include:



- **Biological Nanorobots:** Utilizing biological materials, such as DNA or proteins, to construct nanorobots. These nanorobots may mimic natural biological entities, providing compatibility with biological systems and potential for integration with living organisms<sup>13</sup>.
- **Mechanical Nanorobots:** Constructing nanorobots using traditional materials with mechanical properties at the nanoscale. These robots may have moving parts, gears, or other mechanical components, allowing for precise movements and tasks<sup>13</sup>.
- **Hybrid Nanorobots:** Combining biological and synthetic components to create hybrid nanorobots. This approach harnesses the strengths of both biological and mechanical systems, offering enhanced functionality and adaptability<sup>14</sup>.
- **Chemical Nanorobots:** Employing chemical processes and reactions at the nanoscale to achieve specific tasks. These nanorobots often rely on molecular interactions for propulsion, sensing, and other functionalities<sup>15</sup>.
- **Swarm Nanorobots:** Deploying a large number of simple nanorobots that work collaboratively to accomplish complex tasks. This approach draws inspiration from swarm intelligence observed in nature, where collective behavior leads to emergent capabilities<sup>14</sup>.
- **Teleoperated Nanorobots:** Controlling nanorobots remotely using external forces, magnetic fields, or other means. This approach allows for precise manipulation and navigation within the body or other environments<sup>15</sup>.
- **Bacteria-Based Nanorobots:** Leveraging living bacteria as the basis for nanorobots, genetically modifying them to perform specific tasks. Bacteria-based nanorobots can be designed for drug delivery, sensing, or other applications<sup>14</sup>.
- **Diamondoid Nanorobots:** Building nanorobots from diamondoid materials, which are nanoscale structures resembling diamond. Diamondoid nanorobots offer exceptional strength and stability, making them suitable for various applications.
- **Biohybrid Nanorobots:** Integrating biological components, such as cells or tissues, with synthetic materials to create biohybrid nanorobots. These entities can leverage the advantages of living systems while incorporating engineered functionalities<sup>15</sup>.
- **Magnetic Nanorobots:** Utilizing magnetic fields for the manipulation and navigation of nanorobots. Magnetic nanorobots can be guided to specific locations within the body or other environments using external magnetic forces<sup>16</sup>.

#### **Ideal Characteristics of Nanorobots:**

They lie in the size range of 0.5 to 3 microns with 1-100 nm parts; else they can block capillary flow.

Nanorobots prevent itself from being attacked by the immune system by having a passive, diamond exterior.

It has the capability to communicate with the doctor by encoding messages to acoustic signals at carrier wave frequencies of 1-100 MHz.

It can produce multiple copies of it to replace worn-out units, a process called self-replication<sup>17</sup>.

#### **COMPONENTS OF NANOROBOTS**

Various components in nanorobots includes

##### **Nanobearings and nanogears**

Ball-and-stick models or space-filling models can be used to depict bearings and gears. In this instance, silicon, oxygen, and hydrogen are also included in the 206 carbon atoms that make up the bearing. A tiny shaft with a diameter of 2.2 nm is used in the design, and it spins inside of a ring to reduce friction and enable smooth shaft rotation. The molecular gear mechanism is a small, 4.3 nanometer-diameter structure that is part of the nanoscale. It is made with a silicon shell ended with sulfur atoms, has 12 moving components, and a molecular weight of 51009.84 Daltons. These specifics highlight its complexity and possible nanotechnology uses for precise molecular-level processes<sup>18</sup>.



### **Payload**

The nanorobot carries a drug in a void section and, when inside the body, releases the drug at the targeted site of action<sup>18</sup>.

### **Micro-camera**

Nanorobots move through the circulatory system by the control and monitoring system, which also keeps checks on their performance.

### **Electrodes**

Electrodes that can use ions from biological fluids to build a battery can be used for both tumor therapy and other purposes. The precise electric currents delivered by protruding electrodes might target and perhaps destroy tumor cells. Surface electrodes, meanwhile, provide regulated application for the elimination of tumors. Although further study and medical oversight are needed for its safety and effectiveness, this novel strategy shows promise<sup>19,20</sup>.

### **Laser**

Hypothetical nanorobots with integrated lasers are utilized in medical procedures including the removal of blood clots and plaque in the arteries. The laser systems on these nanorobots would be able to accurately target and remove blood artery barriers, so restoring blood flow<sup>21</sup>.

### **Ultra-sonic signal generator**

It's a novel idea in the world of medical therapy to employ ultrasonic noises along with nanorobots to remove kidney stones. With this strategy, nanorobots that can travel the urinary system and locate kidney stones would be used. The use of ultrasonic waves might help the body naturally eliminate kidney stones by dissolving or breaking them down into smaller, more readily accessible bits<sup>20</sup>.

### **Swimming tail**

Nanorobots use a variety of transportation techniques to penetrate the body. They move using motors and, occasionally, manipulator arms or mechanical legs for greater mobility in the challenging environment. Their control systems are run by specialist software created to mimic and govern the behavior of nanorobots in fluid environments with Brownian motion's unpredictable motions. Nanorobots are extremely useful for medical applications since they can identify target molecules thanks to their chemical sensors. Ant colony optimization (ACO), artificial bee colony (ABC), and particle swarm optimization (PSO), three main swarm intelligence techniques, are harnessed to enhance their collective intelligence and adaptability, promising to revolutionize healthcare by performing precise and effective tasks inside the human body. Nevertheless, as this technology develops, addressing safety, control, and ethical issues becomes essential<sup>21,22</sup>.

### **Nanocomputers**

In fact, having some kind of internal processors is essential for effective nanorobot activity within the human body. These tiny computer systems act as the brains of nanorobots, giving them the ability to carry out precise, focused activities inside the complex atmosphere of the human body. These onboard computers provide control and real-time monitoring of these small devices, making them useful tools for doctors. Physicians can improve the accuracy and effectiveness of medical treatments by programming nanorobots to go through blood veins or tissues, deliver medications to precise areas, or repair damaged cells. These computers' connectivity capabilities allow medical professionals to monitor the progress of the nanorobots and make necessary changes<sup>23</sup>.

### **Application**

- **Drug Delivery:** Nanorobots are envisioned to revolutionize drug delivery by providing precise and targeted transport of therapeutic agents within the body. With their nanoscale size, these robots can navigate through the bloodstream, avoiding obstacles and delivering medications directly to specific cells or tissues. This targeted drug delivery approach reduces systemic side effects and enhances the therapeutic impact of drugs, particularly in treating diseases like cancer where localized treatment is critical<sup>24</sup>.
- **Diagnostics:** Nanorobots equipped with advanced sensors and imaging capabilities have the potential to transform diagnostics. They can be designed to navigate through the body, detecting molecular biomarkers



associated with various diseases. By providing real-time, high-resolution data, nanorobots offer the opportunity for early and accurate disease diagnosis, paving the way for timely interventions and personalized treatment strategies<sup>25</sup>.

- **Cancer Treatment:** Nanorobots hold significant promise in cancer treatment by overcoming the challenges associated with conventional therapies. These robots can deliver anticancer drugs directly to tumor sites, enhancing treatment efficacy while minimizing damage to healthy tissues. Additionally, nanorobots can be engineered to identify cancer cells and facilitate targeted therapies, making them invaluable tools in the fight against cancer<sup>26</sup>.
- **Minimally Invasive Surgery:** Nanorobots enable minimally invasive surgical procedures by navigating through the body with precision. These robots can perform intricate tasks at the cellular or molecular level, reducing the need for large incisions and promoting faster recovery times. Their ability to access hard-to-reach areas and execute precise maneuvers makes them instrumental in advancing surgical techniques.
- **Regenerative Medicine:** Nanorobots contribute to regenerative medicine by participating in tissue repair and regeneration processes. They can be designed to deliver growth factors or manipulate cellular activities to enhance healing and regeneration. This application has implications for treating injuries, degenerative diseases, and promoting overall tissue health<sup>27,28</sup>.
- **Continuous Monitoring and Intervention:** The continuous monitoring capabilities of nanorobots allow for real-time assessment of physiological parameters within the body. These robots can provide valuable data on health metrics, enabling proactive healthcare management. Moreover, nanorobots can intervene autonomously or under external guidance to address emerging health issues, such as blood clot removal or targeted tissue repair<sup>29</sup>.
- **Environmental Cleanup:** Nanorobots can be deployed in environmental cleanup efforts to address pollution and remove contaminants at the molecular level. Their ability to selectively target and neutralize pollutants offers a promising solution for remediating environmental damage, contributing to sustainability efforts<sup>27,30</sup>.
- **Nano-Scale Construction:** In the realm of nanotechnology and materials science, nanorobots play a role in assembling nanoscale structures with precision. This application contributes to the development of advanced materials, such as nanocomposites and nanoelectronics, with tailored properties and functionalities<sup>25,31</sup>.
- **Targeted Imaging:** Nanorobots equipped with advanced imaging technologies provide detailed and real-time imaging of biological structures. This application is particularly valuable in medical diagnostics and research, offering insights into cellular and molecular processes with unprecedented clarity<sup>32</sup>.
- **Neurological Applications:** Nanorobots can be applied in neurological contexts for targeted drug delivery or monitoring within the brain. Their ability to cross the blood-brain barrier opens avenues for treating neurological disorders and gaining insights into brain function at the cellular level<sup>33</sup>.
- **Blood Clot Removal:** Navigating through the bloodstream, nanorobots can target and remove blood clots or blockages. This application is crucial for preventing and treating cardiovascular issues, offering a minimally invasive alternative to conventional interventions<sup>25</sup>.
- **Vaccine Delivery:** Nanorobots contribute to vaccine delivery by enhancing the stability and targeted release of vaccine components. This precision in vaccine administration can lead to more effective immunization strategies, especially in the development of vaccines for emerging infectious diseases<sup>26</sup>.
- **Agriculture and Pest Control:** Nanorobots find applications in agriculture for targeted delivery of pesticides or fertilizers at the cellular level in crops. This targeted approach minimizes environmental impact, enhances agricultural efficiency, and reduces the overall use of agrochemicals.
- **Water Purification:** Nanorobots can assist in water purification by selectively targeting and removing contaminants, ensuring the production of clean and safe drinking water. This application has the potential to address challenges related to water pollution and improve water quality worldwide.
- **Nanorobots in heart surgery:** The integration of nanorobots in heart surgery represents a paradigm shift in cardiovascular medicine, offering unprecedented opportunities for precision, minimally invasive interventions,





and enhanced patient outcomes. One notable application lies in precision drug delivery, where nanorobots can navigate through the intricate network of blood vessels to deliver medications directly to targeted areas within the heart<sup>35</sup>. This approach minimizes systemic side effects and improves the efficacy of treatments for conditions like ischemia or arrhythmias. Nanorobots also hold immense potential in addressing arterial blockages by navigating through blood vessels and removing plaque, providing a less invasive alternative to traditional interventions such as angioplasty. With advanced imaging capabilities, nanorobots enable real-time visualization of the heart's internal structures, aiding surgeons in precisely navigating and performing procedures. In cases of blood clot formation within the heart, nanorobots could be designed to target and remove these clots promptly, mitigating the risk of complications. Furthermore, nanorobots may play a pivotal role in facilitating tissue repair and regeneration in the aftermath of cardiac injuries, delivering growth factors or stem cells to specific regions<sup>34</sup>. Continuous intracardiac monitoring by nanorobots enhances the understanding of dynamic physiological parameters, allowing for proactive management of cardiac conditions. As research progresses, the incorporation of nanorobots into heart surgery holds the promise of ushering in a new era of personalized, targeted, and minimally invasive interventions, ultimately improving the overall landscape of cardiovascular care<sup>36</sup>.

#### **Application of Nanorobots in Dentistry:**

- **Precision Drug Delivery:** Nanorobots can be employed for targeted drug delivery in dentistry. By navigating through the oral environment, nanorobots can deliver therapeutic agents precisely to specific dental tissues, such as the periodontal ligament or dental pulp, for localized treatment of infections or inflammation<sup>37</sup>.
- **Oral Disease Detection and Diagnosis:** Nanorobots equipped with advanced sensors can assist in the early detection and diagnosis of oral diseases. They can navigate through the oral cavity, identifying biomarkers associated with conditions like periodontal disease or oral cancers, providing dentists with real-time diagnostic information.
- **Biofilm Removal:** Nanorobots can contribute to more effective biofilm removal on teeth surfaces. These robots can target and disrupt bacterial biofilms, preventing the formation of plaque and reducing the risk of dental caries and periodontal disease.
- **Tooth Repair and Regeneration:** Nanorobots may play a role in tooth repair and regeneration by delivering growth factors or bioactive molecules to damaged dental tissues. This could potentially enhance the natural healing process and promote the regeneration of dentin or enamel<sup>38</sup>.
- **Orthodontic Applications:** In orthodontics, nanorobots could assist in precise tooth movement by delivering therapeutic agents or exerting controlled forces at the cellular level. This could contribute to more efficient and comfortable orthodontic treatments.
- **Surgical Procedures:** Nanorobots can be employed in minimally invasive surgical procedures within the oral cavity. For example, they could assist in precise tissue removal during dental surgeries or facilitate the placement of dental implants with enhanced accuracy.
- **Dental Imaging Enhancement:** Nanorobots equipped with imaging capabilities can provide detailed imaging of dental structures. This could improve diagnostic procedures, allowing for more accurate assessments of dental conditions and treatment planning<sup>39</sup>.
- **Sensory Applications:** Nanorobots may be designed with sensory capabilities to monitor changes in the oral environment. This includes detecting pH levels, microbial activity, or early signs of dental erosion, enabling preventive interventions and personalized oral care.
- **Salivary Diagnostics:** Nanorobots could be utilized for salivary diagnostics, analyzing saliva for biomarkers associated with various systemic and oral health conditions. This non-invasive approach could offer valuable insights into a patient's overall health.



- **Customized Dental Materials:** Nanorobots may contribute to the development of customized dental materials with improved properties. This includes the creation of nanocomposites for dental fillings or coatings that offer enhanced durability and antimicrobial properties<sup>40,41</sup>.

#### **Future Scope**

##### **Targeted Drug Delivery:**

- Advancements for more precise and efficient drug delivery.
- Enhanced payload capacity and adaptive release mechanisms.

##### **Diagnostics and Imaging:**

- Improved imaging capabilities for early and accurate disease detection.
- Molecular and cellular-level diagnostics using nanorobots.

##### **Minimally Invasive Surgery:**

- Evolution of surgical techniques with greater precision.
- Faster recovery and reduced disruption to healthy tissues.

##### **Regenerative Medicine:**

- Targeted delivery of regenerative therapies for tissue repair.
- Applications in treating injuries, degenerative diseases, and organ damage.

##### **Real-Time Monitoring and Intervention:**

- Continuous monitoring of physiological parameters.
- Early detection for timely interventions and personalized healthcare.

##### **Neurological Applications:**

- Crossing the blood-brain barrier for targeted neurological treatments.
- Potential applications in neurodegenerative diseases and brain tumors.

##### **Immunotherapy Enhancement:**

- Precise delivery of immune-modulating agents for improved response.
- Advancements in cancer immunotherapy.

##### **Integration with IoT and AI:**

- Interconnected systems for remote monitoring and data analysis.
- Automated responses and more efficient healthcare delivery.

##### **Biohybrid Systems:**

- Development of nanorobots combining synthetic and biological components.
- Synergizing advantages of both systems for innovative medical applications.

##### **Clinical Trials and Translation:**

- Extensive clinical trials to validate safety and efficacy.
- Translation of nanorobot technologies from lab to clinical practice.

##### **Addressing Challenges:**

- Overcoming biocompatibility and safety concerns.
- Ethical considerations and regulatory framework development.

##### **Interdisciplinary Collaboration:**

- Continued collaboration between nanotechnology, medicine, and engineering.
- Ensuring responsible and collaborative development in healthcare.

## **II. CONCLUSION**

The applications of nanorobots in medical techniques hold immense potential for transformative advancements. From targeted drug delivery and precise diagnostics to minimally invasive surgery and regenerative medicine, nanorobots offer promising solutions for personalized and efficient healthcare. Real-time monitoring, enhanced immunotherapy,



and integration with IoT and AI further broaden their impact. Addressing challenges and fostering interdisciplinary collaboration are crucial for responsible development, ensuring a future where nanorobots play a pivotal role in advancing medical practices.

### REFERENCES

- [1]. Robert Freitas Jr A. The Ideal Gene Delivery Vector: Chromalloytes, Cell Repair Nanorobots for Chromosome Replacement Therapy, *Journal of Evolution and Technology*, 16(1), 2007, 1-97.
- [2]. Pingoud A, Jeltsch A. Structure and function of type II restriction endonucleases, *Nucleic Acids Res*, 29(18), 2001, 3705-3727.
- [3]. Kraulis P J, Raine A R, Gadhavi P L, Laue E D. Structure of the DNA-binding domain of zinc GAL4, *Nature*, 356(6368), 1992, 448-450.
- [4]. Von Hippel P H, Berg O G. Facilitated target location in biological systems, *J. Biol. Chem*, 264(2), 1989, 675-678.
- [5]. Robinson C R, Sligar S G. Changes in solvation during DNA binding and cleavage are critical to altered specificity of the EcoRI endonuclease, *Proc. Natl Acad. Sci, (USA)*, 95(5), 1998, 2186-2191.
- [6]. Sidorova N Y, Rau D C. Differences in water release for the binding of EcoRI to specific and nonspecific DNA sequences, *Proc. Natl Acad. Sci, (USA)*, 93(22), 1996, 12272-12277.
- [7]. Jeltsch A, Alves J, Wolfes H, Maass G, Pingoud A. Pausing of the restriction endonuclease EcoRI during linear diffusion on DNA, *Biochemistry*, 33(34), 1999, 10215-10219.
- [8]. Grundt K, Skjeldal L, Anthonsen H W, Skauge T, Huitfeldt H S, Ostvold A C. A putative DNA-binding domain in the NUCKS protein, *Arch. Biochem. Biophys*, 407(2), 2002, 168-175.
- [9]. Shioi S, Ose T, Maenaka K, Shiroishi M, Abe Y, Kohda D, Katayama T, Ueda T. Crystal structure of a biologically functional form of PriB from *Escherichia coli* reveals a potential single -stranded DNA-binding site, *Biochem. Biophys. Res. Commun*, 326(4), 2005, 766-776.
- [10]. Trumbore M W, Wang R H, Enkemann S A, Berger S L. Prothymosin alpha in vivo contains phosphorylated glutamic acid residues, *J. Biol. Chem*, 272(42), 1997, 26394-26404.
- [11]. Dobbstein M, Strano S, Roth J, Blandino G. p73-induced apoptosis: a question of compartments and cooperation, *Biochem. Biophys. Res. Commun*, 331(3), 2005, 688- 693.
- [12]. Chapter 24. Mechanics of Biological Nanotechnology, in Bharat Bhushan, ed., *Springer Handbook of Nanotechnology*, Springer-Verlag, Berlin, 2004, 755-760.
- [13]. Jinesh K B, Frenken J W M. Capillary condensation in atomic scale friction: how water acts like a glue, *Phys. Rev. Lett*, 96(16), 2006, 166103.
- [14]. Egger M, Heyn S P, Gaub H E. Synthetic lipid-anchored receptors based on the binding site of a monoclonal antibody, *Biochim. Biophys. Acta*, 1104(1), 1992, 45- 54.
- [15]. Marius P, Alvis S J, East J M, Lee A G. The interfacial lipid binding site on the potassium channel KcsA is specific for anionic phospholipids, *Biophys. J*, 89(6), 2005, 4081-4089.
- [16]. Hornig C, Albert D, Fischer L, Hornig M, Radmark O, Steinhilber D, Werz O. 1 Oleoyl-2-acetyl-glycerol stimulates 5- lipoxygenase activity via a putative (phospho) lipid binding site within the Nterminal C2-like domain, *J. Biol. Chem*, 280(29), 2005, 26913-26921.
- [17]. Van Loon D, Berkhout T A, Demel R A, Wirtz K W. The lipid binding site of the phosphatidylcholine transfer protein from bovine liver, *Chem. Phys. Lipids*, 38(1-2), 1985, 29-39.
- [18]. Lerche M H, Kragelund B B, Bech L M, Poulsen F M. Barley lipid-transfer protein complexed with palmitoyl CoA: the structure reveals a hydrophobic binding site that can expand to fit both large and small lipid-like ligands, *Structure*, 5(2), 1997, 291-306.
- [19]. Springer T A, Wang J H. The threedimensional structure of integrins and their ligands, and conformational regulation of cell adhesion, *Adv. Protein Chem*, 68, 2004, 29-63.
- [20]. Faury G. Role of the elastin-laminin receptor in the cardiovascular system, *Pathol. Biol. (Paris)*, 46(7), 1998, 517-526.





- [21]. De Souza S J, Brentani R. Collagen binding site in collagenase can be determined using the concept of sense-antisense peptide interactions, *J. Biol. Chem*, 267(19), 1992, 13763-13767.
- [22]. Steffensen B, Wallon U M, Overall C M. Extracellular matrix binding properties of recombinant fibronectin type II-like modules of human 72-kDa gelatinase/type IV collagenase. High affinity binding to native type I collagen but not native type IV collagen, *J. Biol. Chem*, 270(19), 1995, 11555-11566.
- [23]. Sumper M, Hallmann A. Biochemistry of the extracellular matrix of *Volvox*, *Int. Rev. Cytol*, 180, 1998.
- [24]. Freitas Jr R A. Exploratory design in medical nanotechnology: A mechanical artificial red cell, *Artificial Cells, Blood Substitutes, and Immobil. Biotech*, 26(4), 1998, 411-430.
- [25]. Prausnitz M R. Microneedles for transdermal drug delivery, *Adv. Drug Delivery Rev*, 56(5), 2004, 581-587.
- [26]. Christoforidis J B, Warden S M, Demattia M, D'Amico D J. In vivo redirection of retinal blood flow into borosilicate microcannulas in rabbits, *Graefes Arch. Clin. Exp. Ophthalmol*, 244(8), 2006, 1022-1025.
- [27]. Olsen T W, Feng X, Wabner K, Conston S R, Sierra D H, Folden D V, Smith M E, Cameron J D. Cannulation of the suprachoroidal space: a novel drug delivery methodology to the posterior segment, *Am. J. Ophthalmol*, 142(5), 2006, 777-787.
- [28]. Li J, Li X, Luo T, Wang R, Liu C, Chen S, Li D, Yue J, Cheng Sh and Sun D: Development of a magnetic microrobot for carrying and delivering targeted cells. *Science Robotics* 2018; 3: 8829.
- [29]. Mehra P and Nabhi K: A Nanorobotics the changing face of dentistry. *International Journal of Science and Research* 2016; 5(3): 192-97.
- [30]. Padovani GC, Feitosa VP, Sauro S, Tay FR, Duran G, Paula AJ and Duran N: Advances in dental materials through nanotechnology facts, perspectives and toxicological aspects. *Trends in Biotechnology* 2015; 33: 621-36.
- [31]. Padovani GC, Feitosa VP, Sauro S, Tay FR, Duran G, Paula AJ and Duran N: Advances in dental materials through nanotechnology facts, perspectives and toxicological aspects. *Trends in Biotechnology* 2015; 33: 621-36.
- [32]. Alexander J H, Hafley G, Harrington R A, Peterson E D, Ferguson Jr T B, Lorenz T J, Goyal A, Gibson M, Mack M J, Gennevois D, Califf R M. Efficacy and safety of edifoligide, an E2F transcription factor decoy, for prevention of vein graft failure following coronary artery bypass graft surgery: PREVENT IV: a randomized controlled trial, *Jama*, 294(19), 2005, 2446- 2454.
- [33]. Kulier A, Levin J, Moser R, RumpoldSeitlinger G, Tudor I C, Snyder-Ramos S A, Moehnle P, Mangano D T. Impact of preoperative anemia on outcome in patients undergoing coronary artery bypass graft surgery, *Circulation*, 116(5), 2007, 471-479.
- [34]. Birkmeyer N J, Charlesworth D C, Hernandez F, Leavitt B J, Marrin C A, Morton J R, Olmstead E M, O'Connor G T. Obesity and risk of adverse outcomes associated with coronary artery bypass surgery, *Circulation*, 97(17), 1998, 1689- 1694.
- [35]. Lukeman Hakkim S, Ranjitha A. Replacement of Heart Bypass Surgery Using Nanorobots: A Review, *International Journal of Novel Research in Engineering and Science*, 2(2), 2016, 7-12.
- [36]. Sharma MK and Gupta R: Nanorobotics the future of medicines. *Research in Pharmacy and Health Sciences* 2016; 2: 51-56.
- [37]. Bi C, Guix M, Johnson B, Jing W and Cappelleri D: Design of microscale magnetic tumbling robots for locomotion in multiple environments and complex terrains. *Micromachines* 2018; 9: 68
- [38]. Mehra P and Nabhi K: A nanorobotics the changing face of dentistry. *The International Journal of Science and Research* 2016; 5(3): 192-97.
- [39]. Cale T S, Lu J Q, Gutmann R J. Threedimensional integration in microelectronics: Motivation, processing, and thermomechanical modelling, *Chemical Engineering Communications*, 195(8), 2008, 847-88.
- [40]. Robert A F J. Current Status of Nanomedicine and Medical Nanorobotics, *Journal of Computational and Theoretical Nanoscience*, 2(1), 2005, 1-25.
- [41]. Kharwade M, Nijhawan M and Modani S. Nanorobots: A Future Medical Device in Diagnosis and Treatment, *Research Journal of Pharmaceutical, Biological and Chemical Sciences*, 4(2), 2013, 1299-1307.

